
S.I.N. 360


DIGITAL WORKFLOW



ACCELERATING INNOVATION

LEAD | SELL | TRAIN | SUPPORT

FULL ARCH DIGITAL WORKFLOW SOLUTIONS





OUR PURPOSE

- LEAD the Industry By Developing the Highest Quality, Most Advanced <u>Technology</u>
- SELL the Most Efficient & Predictable Full Arch Digital Workflow in the Dental Industry
 - TRAIN Surgeons, Dental Office Teams, and Dental Labs on the Latest Workflow for Adoption of the <u>Simplified</u>, Scalable S.I.N. 360 Digital Workflow
- SUPPORT Our Customers with World-Class Care to Ensure the Success and Longevity of Your Practices

FULL ARCH DIGITAL WORKFLOW

EDU & VIDEO DX & DESIGN

IMPLANTS PRODUCTION CONSUMABLES













We provide first-rate training packages designed with your needs in mind. We also create transformative video content to showcase your skills, office facilities, and support of your patients. Our CAD software and scanners are known for their precision, speed, and ability to improve predictability in your results.

The perfect combination for those seeking superior results with high predictability

We provide options for both same-day and permanent prosthesis development.

Our consumables are available in a vast array of shades and translucencies, all paired with exceptional biocompatibility.



ADVANCED EXOCAD COURSE

- ✓ Develop A Comprehensive Understanding of DentalCAD Design for Full Arch
- ✓ Transform Digital Waxups Into 3D Printed Prosthetics On Site
- Learn To Incorporate Photogrammetry Data To Create A Full Arch Design

VIDEO CONTENT

- ✓ Patient Transformation Testimonial & In-Office Tour Video
- 4K Ultra High Definition Video Production

WORKFLOW IMPLEMENTATION

- ✓ Capturing Quality Digital Records & Photos
- ✓ 3D Printing
- ✓ Intra-Oral Scanning for Full Arch Diagnostics
- ✓ Stain & Glaze

Digital Photogrammetry Impressions








SOFTWARE

exocad



EXOPLAN

- Prosthetic Driven Approach to Implant Placement
- Surgical Guide Designs



EXOCAD - DENTALCAD DB

- Screw-Retained Bridges
- Night Guard Design
- Easy Bite Adaptation

- ✓ Immediate Denture Design
- ✓ Digital Wax Up Pre-Planning
- ✓ Smile Creator



DIGITAL SCANNERS

SHINING 3D°



, OC CONTROL OF THE PARTY OF TH

WIRELESS 3D FACE SCANNER

- Fast Scan Speed
- Automatic Alignment of Data
- Seamless exocad Integration
- Easy Attachment Of Photogrammetry Files
- ✓ Seamless Exocad Integration

WIRELESS INTRA-ORAL SCANNER

- Motion Sensing
- ✓ Long-Lasting Battery Life
- Safe and Efficient Scan Experience
- ✓ Highly Efficient Digital Chairside Solution
- Minimal Contact Motion Sensing



MÎCRONMAPPER



MICRONMAPPER - PHOTOGRAMMETRY CAMERA

ELIMINATES:

- Verification Jigs
- Temp Cylinders

Analogs

- Impression Material
- Impression Copings

FEATURES:

- ✓ Intuitive Software
- ✓ Ten Scan Bodies
- Tissue Mapping
 Functionality

- ✓ 400 Uses/ Scanbody
- Instant STL File
- Lightweight ForSingle Handed Use



BONE FIXATION KIT




EACH KIT CONTAINS:

- ✓ Five Gold Bone Screws: 2.0 mm x 10 mm
- Five Gold Bone Screws: 2.0 mm x 12 mm
- One Organization Block
- One Screwdriver & Handle
- One Pilot Drill: PLD 1.4
- 🕢 Two Drivers: IHA 13 (Long & Short)



GOLD: MS-20-010 & 012







PRE-MILLED ABUTMENTS



VERSE PORTFOLIO OF PRE-MILLED BLANKS

One Connection for All Implant Diameters

AVAILABLE FOR MULTIPLE BLANK HOLDERS

- Medentika Holder
- DESS® Holder

MULTIPLE SIZE OPTIONS

- 11.5 MM Diamater
- 16.5 MM Diameter









16° MORSE TAPER CONNNECTION WITH PLATFORM SWITCH

One Connection for All Implant Diameters

HYDROPHILIC NANO SURFACE

- Extensive Research Ensures High Predictability
- Double Acid-Etching of HA Nano Surface Accelerates Osseointegration

AGGRESSIVE DUAL THREAD


Ensures Primary Stability for Same-Day Immediate Loading





ABUTMENT ANGLES:

IMPLANTS



PRH 30 SCREW - 4 PACK



- Anodized Pink
- ✓ 1.4 mm Screw Length
- Pressure Fits on Driver for Secure Connection
- Designed for Direct to MUA Restorations
- ✓ Three Lengths of 1.2 mm Prosthetic Drivers

LAB SUPPORT

We Enable Labs to Design & Manufacture Dental Restorations with Digital Design Courses, Routinely Updated Digital Libraries & CAD/CAM Solutions Built with Precision.



IMMEDIATE PRODUCTION

SprintRay





3D PRINTER SYSTEM

- 20 Minute Print Time
- Resin Approved for Final Restoration
- Pro S Arch Kit
- ✓ High Volume Hub & Spoke Operations

- Enables True Chairside Workflows
- ✓ Fully-Automated, Two-Stage Wash & Dry System.
- OnX Tough 2 available in Bleach, A1 & B1
- Highly-Accurate Model Materials



FINAL PRODUCTION

imes-icore®





MILLING MACHINE

- Integrated Water Supply, Filtration & Compressed Air Supply
- Milling Solutions for Single Unit & Full Arch Cases
- Two Options for C-Clamps
- ✓ High-Performance Spindle

C-CLAMPS

- Enables 5-Axis Machining up to 90 Degrees
 Axis Adjustments
- Front Side Machining
- Extension of the Technical Possibilities of your Milling Machine
- Ideal for Practices Specializing in Full-Mouth Rehabilitation






SINTERING FURNACE

- Best-in-Class Oven Technology
- Remote Opperation Functionality
- ✓ Complete Customization for any Strategy



PORCELAIN FURNACE

- Excellent Firing Results Across All Temperature Ranges
- Remote Opperation Functionality
- Auto-Dry Function for Precise Drying as well as Very Fast and Targeted Heating and Cooling



CONSUMABLES



TITANIUM DISCS

✓ Grade Five & Two





ONX TOUGH RESIN

- Hybrid Ceramic Resin Produces Highly

 Aesthetic Results
- Eliminates the Need for Any Mixing or Rolling



PEEK & PMMA DISCS

Shaded & Multicolor Available



ZIRCONIA PUCKS

- Highly Compressed Material Ensures Precise

 Margins & Edges
- ✓ Fine Layering & Homogeneous Coloring
- ✓ Consistently Reproducible Shades

Epikut^s

EPIKUT S PLUS WAS IDEALIZED FOR THOSE WHO WANT TO REDEFINE THE CONCEPT OF DENTAL IMPLANTS.


With a cutting and compressive design, double inverted support screws, combined with the ultra-thin surface Plus, which is produced by double acid-etching followed by application of a hydroxyapatite coating HAnano.


Learn More





TIMELINE





S.I.N. 360 SELL TRAIN SUPPORT

VISIT OUR WEBSITE



WWW.SIN360.US | 866.507.9315